

## New Willi Fox rapid test

# COVID-19 IgM/ IgG rapid test

From day 7 or day 8 IgM and IgG antibodies to SARS-CoV-2 are detectable!

With the new Willi Fox COVID-19 IgM / IgG you can find out, whether a SARS-CoV-2 infection has already occurred and whether there is already an immunization.



### Detection is possible with:

- Whole blood (e.g.finger blood)
- Serum
- Plasma

### Performance IgM:

Relative sensitivity:	93,7% (86,0%-97,3%)*
Relative specificity:	99,1% (96,8%-99,8%)*
Overall agreement:	97,7% (95,4%-98,9%)*

\*95% Confidence Interval

### Performance IgG:

Relative sensitivity:	98,8% (93,5%-99,8%)*
Relative specificity:	98,7% (96,2%-99,5%)*
Overall agreement:	98,7% (96,7%-99,5%)*

\*95% Confidence Interval

## Quick and reliable

Long shelf life and easy storage!